



Thawing and warming of plasma and blood products – dry, safe and reliable

Barkey plasmatherm

The Barkey plasmatherm thaws Fresh Frozen Plasma (FFP), Plasma Frozen within 24 hours after Phlebotomy (PF24), Cryoprecipitated AHF, Whole Blood, Red Blood Cells (RBC) and crystalloid infusion solutions.

Key features

- Thaws 4 bags at a time
- Simple and intuitive application
- Dry thawing device = reduced risk of contamination (BSH* recommendation)
- High standard of hygiene due to dry warming
- Smooth surfaces and rounded edges
- Easy to clean through wipe disinfection
- 2 leak detection sensors
- Shortcuts for blood, plasma and user program
- Large display unit
- Replace water once a year
- Warming cushions tightly conform to product

User 10 Seriel no	Sample 10 1 Sample 10 2 Sample 10 2 Sample 10 2 Sample 10 2 Sample 10 4	Sangle ID S Sangle ID 6 Sangle ID 7 Sangle ID 8		IP Address Directory	192.162.100.96 Connect Disconnect
hugun	Sart Sdorn gi Erd Max	Time (bhowman)	Tenpenture		
Event	Description	Temperature	Time (bhommess)		
Connected for					Vector 10

Document the thawing process with the Barkey TCP Logging Tool

^{*} recommended from the Bristish Society for Hematology (BSH): British Society of Haematology Guidelines on the spectrum of fresh frozen plasma and cryoprecipitate products: their handling and use in various patient groups in the absence of major bleeding, in: British Journal of Hematology. 2018, 181, 54 67, p. 69



Barkey plasmatherm

Accessories (optional)

- Documentation software Barkey TCP Logging Tool (via Ethernet connection)
- Barcode reader
- Log printer

Hygiene

- Reduced risk of contamination due to dry thawing and warming (no direct contact with water)
- Easy to clean through wipe disinfection
- Corrosion-resistant

Safety

- 2 redundant overtemperature alarms
- Adjustable temperature range from +30°C to +37°C
- Visual and audible alarm after thawing cycle is complete
- Visual inspection due to clear cushions and lid
- 2 leak detection sensors





Product	Barkey plasmatherm
Power supply	100-120 V 50/60 Hz
Set temperature	+30°C to max. +37°C
Dimensions	13.5 x 23.6 x 12.6 in (WxDxH)
Weight without filling	42 lbs
Filling quantity	9 liters

MANUFACTURED BY:

Barkey GmbH & Co. KG Technical amendments and errors excepted.

Interfaces/Devices	Barcode scanner, Log printer,
(optional)	Ethernet port
MDR classification	lla
Safety class	I, protective earthing
Certification	CE 0123, FDA, UL
Barkey is certified	DIN EN ISO 13485

DISTRIBUTED BY:

Barkey Corporation 500 West Cummings Park Suite 1900 Woburn, MA 01801 1 866 BARKEY3 info@barkey-us.com www.barkey-us.com