

From Azenta Life Sciences



Barkey S-line C&G

Reliable Fluid Path Warmer

- Warms fluids like cell culture media as well as blood and cellular products used in laboratories and manufacturing.
- Alerts user with a visual and acoustic alarm in the event of an error.
- One preset temperature for easy operation.
- One-touch control and no additional disposables required.



Safety

- Integrated temperate sensors
- Overtemperature alarms
- Smooth surfaces allow for easy disinfecting procedures
- Provides operational safety feature due to preset target temperature value

Application

A safe and reliable system which helps support the warming of cellular products through all steps of tech transfer. The Barkey S-line C&G is applicable for C&G Research & Development up to the point of manufacturing.

Wide range of services

- Installation Qualification (IQ)
- Operational Qualification (OQ)
- Performance Qualification (PQ)
- Preventive Maintenance (PM)

MANUFACTURED BY:

Barkey GmbH & Co. KG Technical amendments and errors excepted.



DISTRIBUTED BY:

Barkey Corporation 500 West Cummings Park Suite 1900 Woburn, MA 01801 1 866 BARKEY3 info@barkey-us.com www.barkey-us.com