

pi-ZeZ-EW-UUU /-U3 viaiguard pro Flyer U3

Key features

- Barkey vialholder ensures proper thawing process (customizable for multiple types of vials)
- Ensures maximal product recovery for small volume thawing
- Processing multiple or single volumes at a time
- Vertical thawing as preferred thawing method
- GMP and GLP compliant
- Optimal for QC needs
- Thightly conforms to most commercially available cryovials and vials
- Versatile and low maintenance device
- For reproducible standardized process, thaw after thaw

Reproducibility

- Automated and robust process
- Customized settings
- Reduced risk of human error

Services

- Qualification package IQ, OQ, PQ
- Maintenance plan PM
- Expert training for users

Accessories (optional)

- Documentation software Barkey TCP Logging Tool (via Ethernet connection)
- Bar code reader



vialguard pro

⇒Barkey

0

MANUFACTURED BY:

Barkey GmbH & Co. KG Technical amendments and errors excepted.



DISTRIBUTED BY:

→ Barkey vialguard pro

Barkey Corporation 500 West Cummings Park Suite 1900 Woburn, MA 01801 1 866 BARKEY3 info@barkey-us.com www.barkey-us.com